
Windows Task scheduler for ELRImporter

On Windows, Task Scheduler is a tool that allows us to create and run tasks automatically.


To create a task, use the following steps:

- · Open Start.
- Search for **Task Scheduler**, and click the top result to open the app.




Click the **Action menu**.

Select the "Create Task" option.




General tab:


In the "Name" field, type a short descriptive name for the task – for example, ELRImporterTask. Select the "Run whether user is logged on or not" option.



Triggers tab:

Configure the trigger values for the ELR Importer, as shown in the screen below.





Actions tab:


Use the "Action" drop-down menu and select the **Start a program** option.


Under the "Settings" section, in the "Program/script" field, specify the path for the application.

Program/script: D:\wildfly-10.0.0.Final\nedssdomain\Nedss\BatchFiles\ELRImporter.bat

Add arguments(optional): >ElrImporter.output


Start in(optional): D:\wildfly-10.0.0.Final\nedssdomain\Nedss\BatchFiles






Conditions tab:

Use this option with the default selections.



Settings tab:

Use this option with the default selections.



History tab:

Under the **Actions** section on the right side, there is an option to Enable/Disable All Tasks History.

